



ClariFlex™

Foldable IOL With OptiEdge

Model CLRFLXB

**ClariFlex™
Silicone Biconvex
Posterior Chamber
CLRFLXB 6 mm
Noninvestigational
Mod C**



*All dimensions are nominal and are specified in mm.
†For a 20.0-diopter lens.

DESCRIPTION

OPTIC CHARACTERISTICS

Diameter:	6 mm
Shape:	Biconvex
Material:	SLM-2 (Silicone) With UV Blocker
Manipulation Holes:	None
Powers:	6.0 to 30.0 D
A Constant:	118.0
Theoretical AC Depth:	4.7 mm
Refractive Index:	1.460 (35°C)

HAPTIC CHARACTERISTICS

Length:	13 mm
Style:	Modified C
Material:	Blue Core Polymethylmethacrylate (PMMA) Monofilament
Angle:	10°
Compressibility:	
11.0 mm:	101 mg (0.10 g)
10.0 mm:	177 mg (0.18 g)
Average Weight:	
Air:	19.3 mg
Water:	2.7 mg

CLRFLXB

COMPRESSIBILITY CHART*



*Data on file at AMO, Inc.

Noninvestigational

**For more information,
call 1.800.366.6554**

ClariFlex™
Foldable IOL With OptiEdge

AMO
ADVANCED MEDICAL OPTICS

© 2002 Advanced Medical Optics, Inc., Santa Ana, CA 92705 USA
® and ™ marks owned by AMO, Inc. www.amo-inc.com